

manicare[®]
salon

Magnifying Pore Vacuum

NEW



* Final packaging

Manicare[®] Salon Magnifying Pore Vacuum

Instantly
CLEARER,
Smoother,
Luminous Skin

For smoother, clearer, more luminous skin instantly, the Manicare[®] magnifying pore vacuum delivers noticeable results in every application. It helps you visibly clear congested skin with ease. Includes a in-built 10xmagnification high definition camera linked to your smart phone via a simple App.



Manicare[®] Salon Magnifying Pore Vacuum

Instantly
CLEARER,
Smoother,
Luminous Skin

Gentle
suction for
sensitive
skin

Helps Clean
Up Dead
Skin During
Extraction

Noticeable
Results in
EVERY
Application



FEATURES



10 X magnifying HD built in camera



HD magnifying camera links to your smart phone, via an App



3 intensity settings: Soft, Normal & Strong



Each kit contains 4 different nozzles to suit your needs.

LARGE CIRCULAR NOZZLE

Suitable for most skin types



OVAL NOZZLE

Most targeted suction for stubborn blackheads/T-zone area



SMALL CIRCULAR NOZZLE

Weaker suction for sensitive skin



MICRODERMABRASION NOZZLE

Helps clear dead skin during extraction

(Microdermabrasion nozzle has a mild abrasive surface which helps to clear away tougher debris including dead skin.)



OPERATION INSTRUCTIONS

Power Display

Battery icon will flash when charging



10 x Magnifying HD Camera

Links to Smartphone camera

Level 3 Intensity – STRONG

Press  three times for the highest intensity

Level 2 Intensity – NORMAL

Press  twice for regular strength

Level 1 Intensity – SOFT

Press  once for gentle suction



Switch on/off

Press for 2-3 seconds to power on/off
Standby mode when screen lights up



USB Charging Cord

Get it on Google Play/Download on the App Store

DOWNLOAD THE APP

Download the Manicare® App on your smart phone (Compatible with Android 5.0 or IOS 12.0 and above).

ACTIVATE THE CAMERA

Switch on the device. On your smart phone, select Settings . Wi-Fi . Manicare® network. Open the Manicare® App and follow the prompts.

QUICK START GUIDE

ACTIVATE | CLEANSE | EXTRACT

1. Activate

Activate the camera through the Manicare® App.

2. Cleanse

Cleanse your skin thoroughly with warm water to relax the pores.

For the best results, use Manicare® magnifying pore vacuum immediately after shower, or a steam bath to (raise the temperature of treated area.) This softens the hardened oil or blackhead within your pores for easier extraction.

3. Extraction

Using your preferred nozzle and intensity, glide the vacuum gently over your skin, keeping the nozzle moving.

Expert Tips: Do not stay on the same spot for more than 3 seconds to prevent bruising your skin.



BEAUTY TIPS

- With all beauty tools, trialling on an insignificant place of the body is necessary. Try using each of the nozzles on the inside of the arm testing the intensity of suction prior to using on the face.
- Like most beauty treatments, it is recommended to avoid use if suffering inflamed skin like eczema, sunburnt or active acne as it can cause further irritation. Do not using this product over active breakouts.
- Try using the Pore Vacuum before bed to allow the skin to rest and avoid redness during the day.
- Before using the pore vacuum, you'll want to gently steam your skin either with a warm, damp towel or in the shower. (Although pores don't technically open, this step helps loosen the sebum.) Then, place the pore vacuum on an area of the skin with congested pores—the nose, for example—starting on the lowest setting.
- Instead of pressing the vacuum down on a particularly clogged area and holding the device there, you'll want to move in fairly light, circular motions, or move the head of the instrument around like a wand in a back and forth or motion.
- Do not stay on the same spot for more than 3 seconds, continuously keeping the tool moving to prevent bruising. If congestion didn't ease in the first time, glide over the same area up to 3 times to help it to ease.
- It's also best not to use the device on your entire face; only use it in areas that look or feel particularly congested.
- After use it is normal for skin to appear lightly flushed in treated area.
- After every application, apply an alcohol-free toner to remove residue, followed by apply a cold compress, e.g. Manicare® Massage Globe to reduce redness before continuing your skincare routine

PRECAUTIONS:

Pore vacuums are generally safe to use, but be sure to use appropriate settings depending on your skin. In other words, switching to the highest setting won't necessarily extract more blackheads.

If you have rosacea, telangiectasias [spider veins], reactive/sensitive skin, or you bruise easily, be sure to use the low setting.

As with any product, if you experience pain, stop using the tool right away.

Do not use on children.

If you are susceptible to one or more of the conditions listed below, please seek advice from a medical professional prior to use:

- Skin allergies
- Sunburn
- Opened wounds
- Itchy skin or skin prone to irritation
- Severe acne or active breakouts

If the congestion hasn't been eased, try it again in a couple of days.

Do not overuse the vacuum as you will risk damaging the skin or potentially breaking a capillary.

CARE TIPS

- Thoroughly wash all the used nozzles with warm water.
- Do not wash the handset. Best to wipe with a damp cloth and leave it to dry after each use.
- Store handset and nozzles in the protective original packaging provided.

TROUBLE SHOOTING WITH APP

Where to download the APP?

You can download the Manicare® APP via Google Play or Apple App store for free.

I can't open the APP, why?

The APP is designed for Android 5.0 or iOS 12.0 and above. Please check if your smartphone is using the latest updated operating system.

How to activate the camera?

1. Switch on the device
2. On your smart phone, select Settings > Wi-Fi > Manicare® network.
3. Open the Manicare® App and follow the prompts.

I have activated the APP but it is not connecting, why?

Check if your WiFi Setting has selected the Manciare® network.